

# CAE BabySIM<sup>®</sup>

Small but powerful, validated physiology



Tiny humans depend on us in big ways. That's why CAE Healthcare skipped no details when developing its smallest member of the patient simulation family. CAE BabySIM is an infant-sized, high-fidelity simulator that is powered with advanced physiology for practicing critical care and trauma care. So real, learners are immersed as they experience seesaw breathing, a palpable fontanel, cooing, and crying. BabySIM also supports intubations, chest compressions, chest tube insertions, and IO insertions to prepare learners for the most critical experiences in trauma care. In BabySIM we deliver critical features, because we know it is critical to be prepared.



# Technical Specifications

## Standard Equipment

BabySIM manikin  
Instructor's workstation  
MUSE software  
2 Patients  
4 Simulated Clinical Experiences (SCEs)  
Power and Communications Unit  
4 SCE development licenses  
Manikin carrying case  
TouchPro wireless emulated patient monitor  
Portable air compressor  
Electronic user guide  
CAE Assurance value plan with customer and technical support, Training for Life™ and option to renew



## Optional Equipment

Tablet PC  
Hands-free cable kit  
Manikin tool kit  
Wall Air kit

## Optional Software

Infant Emergencies Learning Modules  
Infant Nursing Learning Modules

## Manikin

25 ¾ inches (65.4 cm)  
16 ¼ lbs (7.37 kg)

## Electrical

Input: 100-240V, 50/60Hz, 2.3A

## Ambient Temperature Range

Operation: 40°F to 104°F

## Humidity

0% to 90% noncondensing

# Key Features

## Airway

- An anatomically realistic upper airway (oropharynx, nasopharynx and larynx), representing that of a three-to six-month-old patient.
- Intubation: orotracheal, nasotracheal, detection and response to right mainstem
- ET Tube, LMA and Combitube
- Bag-valve-mask ventilation

## Breathing

- Measures the presence or absence of carbon dioxide exhalation synchronized with programmed lung sounds. Lung sounds include normal, biphasic, wheeze bronchovesicular, crackles, prolongation, rhonchi, stridor, wheezing
- Seesaw breathing
- Bilateral and unilateral chest rise and fall
- Spontaneous breathing
- Unilateral chest tube insertion with fluid output
- Unilateral needle decompression detected and generates changes in pulmonary mechanics and gas exchange
- Esophageal intubation is fully supported with gastric distension

## Circulation

- Blood pressure measurement by auscultation and palpation
- Bilateral brachial and femoral pulses
- Automatic pulse deficit occurs with decreased systolic blood pressure
- Bilateral brachial and femoral pulses are synchronized with ECG
- Pulses are independently controllable by the instructor for presence/ absence in the case of trauma to a specific extremity

## Cardiac

Accurately simulates a wide variety of hemodynamic conditions and responses

- Live defibrillation, pacing, and cardioversion
- Defibrillation and cardioversion using live defibrillators, energy is automatically quantified and logged
- Pacing (use of hands-free pads), current is automatically quantified and logged
- 3-lead dynamic ECG display
- Cardiac Sounds include:
  - normal
  - S3, S4, S3 and S4
  - Early Systolic Murmur, Mid Systolic Murmur, Late Systolic Murmur, Pan Systolic Murmur, Late Diastolic Murmur

## CPR

- Correct hand placement, depth, and rate of compressions are reflected in physiological feedback rather than virtual target on instructor's workstation



- Adequate chest compressions result in simulated circulation, cardiac output, central and peripheral blood pressures, and carbon dioxide return

## Vascular Access

- IV access supported via left femoral catheter
- IO site on anterior tibia of right leg

## Neurological

- Software driven tri-state pupils and blinking eyes
- Automatic changes in blinking based on inadequate respiratory and cardiovascular conditions
- Palpable anterior fontanel
- Bulging fontanel capability to simulate increased intracranial pressure

## Urological

- Urine output with and without catheter insertion
- Instructor controlled rate and flow of urine output
- Interchangeable genitalia

## Metabolic System

Metabolic features are physiologically modeled within the software and the results are made available on the instructor workstation

- ABG data displayed corresponds accurately and dynamically to the alveolar concentration of CO<sub>2</sub> and O<sub>2</sub>
- Instructor driven simulated metabolic acidosis and alkalosis

## Pharmacology System

- Pharmacology system models automatically and calculates the pharmacokinetics and pharmacodynamics for six intravenous and inhaled medications
- All patient responses to drugs are automatic, dose dependent and follow appropriate time course

## Trauma

- Secretions from eyes, ears and mouth

## Sounds

- Crying/cooing, breath, heart and bowel sounds
- Pre-recorded sounds and can be looped for continuous replay
- Customized sounds and voices via the provided wireless microphone