

# **J300**

# CPAP NASAL MASK

### **TECHNICAL SPECIFICATIONS**

### Weight J300 Nasal Mask

With headgear (M) 66g Headgear 31g

#### **Material**

Cushion Silicone Headgear Nylon and Spandex

#### Resistance

approx.  $0.22 \text{cmH}_2\text{O}$  @ 50L/min. approx.  $0.95 \text{cmH}_2\text{O}$  @ 100L/min.

### **Dimensions mask main body**

(M-Size)

62mm (L) x 61mm (W) x 60mm (H)

Item number (M-Size): J300-M

# Introducing the NEW J300 Nasal CPAP Mask from Drive DeVilbiss Healthcare.

It's no secret that the mask interface is the primary key to successful CPAP therapy. With this in mind the J300 nasal CPAP mask was developed to overcome the common problems of today when choosing the right nasal mask. This stylish new mask is light weight, discreet, user friendly and offers optimum comfort for the patient.

#### **PRODUCT FEATURES**

- Only 4 mask parts for easy assembly, use and cleaning
- Available as a medium sized mask to fit a wide range of facial contours
- Soft silicone cushion
- Quiet in operation
- Directional exhalation port

The J300 includes an elegant soft material headgear providing ultimate comfort and a stylish appearance. Adjustable fastening straps and quick release clips allows the mask to be removed with one hand. Due to the compact size, ease of use and above all comfort, the J300 is the perfect companion for any CPAP device and its patient.

