



SleepCube Standard

...feel refreshed!

SPECIFICATIONS

Dimensions (H x W x L)

10.8 x 16.5 x 17.8 cm

Weight

1.18 kg

CPAP Pressure Range

3 - 20 cm H₂O

Operating

5° C to 40° C; 0 % to 95 % R.H.
non-condensing

Atmospheric Conditions

Sea level to 9,000 feet

Storage & Transportation

-40° C to + 70° C; 0 % to 95 % R.H.
non-condensing

Sound level

< 26 dB(A)

Maximum Power Consumption

65 W from AC power source
(PAP device only)

Electrical Requirements

100 - 240 VAC, 50/60 Hz

The SleepCube standard is a fixed pressure CPAP designed for the effective treatment of OSA. It is equally suitable for use in the home or the hospital and its robust ergonomic design and small size makes it the ideal travelling companion. The SleepCube Standards' multi-voltage capability means that it can be used virtually anywhere in the world and can also be operated from a 12 V DC power source such as those found in a motor home, recreational vehicle or boat.

The option heated humidifier integrates seamlessly into the system and the SleepCube standard is packed full of all the features needed for a restful and refreshing nights sleep.



The simple LED push button display makes the SleepCube easy to use, even in the dark



Back of the SleepCube - showing filter, patient tubing and download outlet



SleepCube with integrated heated humidifier



The heated humidifier is easy to take apart and clean

Drive DeVilbiss International

DeVilbiss Healthcare GmbH
Kamenzer Straße 3
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40
DE +49 (0) 621 / 1 78 98 0
kontakt@devilbisshc.com



www.drivedevilbiss-int.com

